



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electronic Metronome

**Product Code :**  
LBNY-0004-101000141



**Description :**

It is used as disconnection alarm in music teaching, PE, physics and students' extra-curriculum activities, with supply voltage: 220V±10%, 50HZ?  
Audio output power: 0.5W?  
Output amplitude: 0-4V and tempo: 10-290 time/min.  
Dimension (L\*W\*H CM): 14.5\*14.6\*5.3.

**Technical Specification :**

Electronic Metronome



**Laboratory Instrument India**