



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Metronome

Product Code :
LBNY-0004-101000141



Description :

It is used as disconnection alarm in music teaching, PE, physics and students' extra-curriculum activities, with supply voltage: 220V±10%, 50HZ?
Audio output power: 0.5W?
Output amplitude: 0-4V and tempo: 10-290 time/min.
Dimension (L*W*H CM): 14.5*14.6*5.3.

Technical Specification :

Electronic Metronome



Laboratory Instrument India