



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Simple Timer

**Product Code :**  
LBNY-0004-101000148



**Description :**

It is used to measure time and count up in middle school physics teaching demonstration experiment and is common used with oblique trough orbit.  
It adopts infrared transmitting and receiving device and IC internal device.

**Technical Specification :**

Dimension (L\*W\*H CM): 19.6\*17\*9.4



**Laboratory Instrument India**