



Email : sales@laboratoryinstrumentindia.com

Product Name :
Low Frequency Signal Generator

Product Code :
LBNY-0004-101000154



www.laboratoryinstrumentindia.com

Description :

It is used for electro-physics and electronic technology teaching and research in colleges, secondary vocational schools and various schools, as well as, it is used as special purpose experimental econometric tools in scientific research units. It can simultaneously produce sine wave, triangular wave and square wave, with digital display output signal frequency.

Product dimension (L*W*H CM): 24.5*19*10.8

Technical Specification :

Low Frequency Signal Generator



Laboratory Instrument India