



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Low Frequency Signal Generator

**Product Code :**  
LBNY-0004-101000154



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It is used for electro-physics and electronic technology teaching and research in colleges, secondary vocational schools and various schools, as well as, it is used as special purpose experimental econometric tools in scientific research units. It can simultaneously produce sine wave, triangular wave and square wave, with digital display output signal frequency.

Product dimension (L\*W\*H CM): 24.5\*19\*10.8

**Technical Specification :**

Low Frequency Signal Generator



**Laboratory Instrument India**