




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
GSM and GPRS Telemetry System

**Product Code :**  
LBNY-0002-860003



**Description :**

GSM and GPRS Telemetry System

**Technical Specification :**

GSM / GPRS Telemetry System

Module: Quad band GSM 850/900/1800/1900 MHz

Output Power: SIM Interface: External standard SIM 1.8V & 3V

SMS: Supports Text, FTP and PDU mode

Data: Data circuit Asynchronous and non-transparent up to 14.4 kbps

GPRS Packet Data Feature: Coding Scheme CS1 to CS4, Embedded TCP/IP Stack

Radio Modem: Works with VHF and UHF user adjustable ranges

AT Commands Interface: GSM /GPRS / TCP-IP/FTP

Connectors/Switches: RS-232 Serial or USB, SIM Card Holder, DC power connector  
Antenna connector

Power Supply: Input 6 - 30V DC with reverse connection protection

Enclosure: Water proof and heat resistant material

Mounting: Poll, Wall Mounting

Temperature Operating: - 10°C ~ +70°C

Software: The software in English, compatible with Windows, Installation and operation manuals is included

Antenna: Rubber cover.



**Laboratory Instrument India**