



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Transistor Millivoltmeter

**Product Code :**  
LBNY-0004-101000158



**Description :**

It is a metrical instrument with high precision used in measuring sine voltage. It uses high resistance FET as input resistor converter, and uses multilevel direct coupling amplifier circuit as signal enlarging. It is specially used for electricity experiment in schools physics teaching, and provides professional metrical instrument for factory and scientific research units.

**Technical Specification :**

Dimension (L\*W\*H CM): 26.8\*16\*12.8



**Laboratory Instrument India**