



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Three-Wire Electronic Switch

**Product Code :**  
LBNY-0004-101000160



**Description :**

It can connect oscilloscope and display three signals waveform, be used in school class teaching demonstration and student extracurricular experiment, and be especially suitable for demonstrating three-phase AC waveform in AC curriculum.

**Technical Specification :**

Dimension (L\*W\*H CM): 20\*17\*10



**Laboratory Instrument India**