



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Teaching Q Table

**Product Code :**  
LBNY-0004-101000162



**Description :**

It is used to measure primary and secondary winding of power transformer, power filter choke, speaker crossover, line output transformer of TV Set, field and line deflection yoke, inductance value and Q value of radio cassette bias oscillating coil and current loops.

**Technical Specification :**

Dimension (L\*W\*H CM): 32.3\*16.1\*15.6



**Laboratory Instrument India**