



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Conductor-Insulator Checker

**Product Code :**  
LBNY-0004-101000163



**Description :**

It is used to test conductor and insulator, with acoustic-optic indication. If the green light is switched on and there is a sound, it is a conductor, vice versa, insulator.

**Technical Specification :**

Dimension (L\*W\*H CM): 13.2\*6.6\*8.5



**Laboratory Instrument India**