



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sensor Applied Experimental Apparatus

**Product Code :**  
LBNY-0004-101000180



**Description :**

It is used in thermistor, characteristic experiment, photoresistance characteristic experiment, photoelectric counting basic principle experiment and simple temperature automatic control experiment.

**Technical Specification :**

Dimension (L\*W\*H CM): 18.8\*13.2\*6



**Laboratory Instrument India**