



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spherical Conductors

Product Code :
LBNY-0004-101000192



Description :

The sphere uses copper air-core sphere, with nickel plating on the surface, 94mm in diameter of sphere, 105mm total height between insulating support rod and base, 10.5mm in diameter of support rod, 103mm in bottom diameter of base. It is used to demonstrate static induction and charge distribution of charged conductor. It consists of sphere, insulating support rod and base.

Technical Specification :

Spherical Conductors



Laboratory Instrument India