



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Air Compression Ignition Apparatus

**Product Code :**  
LBNY-0004-101000193



**Description :**

It consists of cylinder, base, end cap, piston and seal ring. The cylinder is made of organic glass, with length of cylinder 130mm, outer diameter 25mm and inner diameter 10mm. diameter of base is 65mm, diameter of handle is 40mm and diameter of piston is 8mm. It is used to demonstrate experiment of "temperature increases with air is compressed".

**Technical Specification :**

Dimension (L\*W\*H CM): 6.5\*23.3



**Laboratory Instrument India**