



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electricity and Electronics Fundamentals Trainer

Product Code :
LBNY-0002-620010



Description :

Electricity and Electronics Fundamentals Trainer

Technical Specification :

Electricity and Electronics Fundamentals Trainer covering the concepts, theory and applications of basic electricity and electronics that are required for technicians.

This comprehensive course is presented in a systematic way starting the student out at the beginning covering basic theory and measurements with meters.

Base Station: Student Base Station: This unit serves as the platform for carrying out the student experiments. Through 2 mm sockets and test leads.

It includes a variety of components such as the following:

Trainer board, component and Circuit based PCB Module

Base Station: Fixed regulated DC

Power Supplies: +5V, +15V, -15V

Variable Regulated DC Power Supply: 0 - 25V

Potentiometers: 4 units

Pulse Switches: 2 units

8 bits LED display.

Toggle Switches: 5 spdt

Speaker: 2.25 inch, 8 ohm/0.25w

One set of 2 mm stackable test leads to perform all required inter-connections.

Set of PCB Board Experiment Modules

Electricity and Electronics Fundamentals Trainer, Electricity and Electronics Fundamentals Trainer Equipments, Electricity and Electronics Fundamentals Trainer Tools, Electricity and Electronics Fundamentals Trainer Tool Kits, Electricity and Electronics Fundamentals Trainer Manufacturers, Electricity and Electronics Fundamentals Trainer Suppliers from India, China, Kenya



Laboratory Instrument India