



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Liquid Internal Pressure Demonstrator

**Product Code :**  
LBNY-0004-101000199



**Description :**

It consists of round box, metal rod, fixed clamp, and upper grooved wheel with knob, rubber belt, connected round box and latex tube with minimal piezometer. It is used to verify that inner pressure of liquid varies with the depth of object; the pressure to every direction is equidistant at the same depth place of liquid inner.

**Technical Specification :**

Liquid Internal Pressure Demonstrator



**Laboratory Instrument India**