



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pointer Electroscopes

**Product Code :**  
LBNY-0004-101000200



**Description :**

It consists of circular base, transparent glass, metal rod, sleeve, pointer, horizontal axis and housing. It can determine the potential difference between two charged bodies, test whether the object is charged, test that the charged object has positive electricity or negative electricity, test the insulating of object and perform the experiment of electrification by influence.

**Technical Specification :**

Dimension (L\*W\*H): 29\*7



**Laboratory Instrument India**