




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wave Synthesis Demonstrator

**Product Code :**  
LBNY-0004-101000201



**Description :**

It uses integrated circuit control and displays waveform with LED. Student can see directly the synthesis of two waves with in-phase, same frequency and different amplitude, as well as, it can demonstrate the synthesis of special form (two reverse waves) of standing wave.

**Technical Specification :**

Dimension (L\*W\*H CM): 50.3\*33\*6.7



**Laboratory Instrument India**