



Email : sales@laboratoryinstrumentindia.com

Product Name :
Resistance Law Demonstrator

Product Code :
LBNY-0004-101000203



Description :

It is used to demonstrate metal conductor law of resistance in middle school physics teaching. It has copper, iron and nickel chrome metal wire. It can qualitatively describe relation between conductor material and resistance, verify quantitative relation that the resistance of a conductor varies directly with its length, conversely with its cross-sectional area.

Technical Specification :

Resistance Law Demonstrator



Laboratory Instrument India