



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Overweight Weightlessness Demonstrator

**Product Code :**  
LBNY-0004-101000232



**Description :**

The dynamometer maximum limit is 2N, the minimum scale is 0.02N. The slot is 3\*50g.  
The outside dimension is 250mm\*160mm\*160mm.

**Technical Specification :**

Overweight Weightlessness Demonstrator



**Laboratory Instrument India**