



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Bed Base

**Product Code :**  
LBNY-0001-820001



#### **Description :**

Psychiatric Bed Base, Psychiatric Bed Base Equipments, Psychiatric Bed Base Tools, Psychiatric Bed Base Tool Kits, Psychiatric Bed Base Manufacturers, Psychiatric Bed Base Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara

#### **Technical Specification :**

7 pressed steel mattress base  
Holes on each mattress base for restraining belt  
Strong bed support (Four legs and additional pipe leg in the middle)  
Base a one piece seamless rotationally moulded design (for easy cleaning and disinfection)  
Constructed from Polyurethane or material of similar characteristic. Bidder to specify Unit fire retardant properties  
Build in mattress retainer to prevent mattress from sliding/shifting  
Base the features to collect fluids (spilled/ during cleaning)  
Static weight bearing of 500kg  
Dimensions: L 2000 x W 1000 x H 500mm +/-200mm  
Weight: 40 kg +/- 10k

---

Base have Omm clearance from ground all-round/no legs/no castors features to allow base to be anchored to the floor (using bolts/screws), fasteners cancelled from patients view.  
Domed/ rounded edges to ensure safety of patients and caregivers.



**Laboratory Instrument India**