



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Shaft Model

**Product Code :**  
LBNY-0004-101000236



**Description :**

It is composed of wheel( including the 120mm's large wheels and 60mm's small wheel), hob, iron shaft, connection blocks, the rocker arm, and the balance block etc.  
It is used to demonstrate the lever principle of the axle.

**Technical Specification :**

Shaft Model



**Laboratory Instrument India**