



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bar Box Dynamometer 1N2N5N10N

**Product Code :**  
LBNY-0004-101000245



**Description :**

It is used to measure the size of the force.  
Specifications: 10 N \* 0.2 N, 5 N \* 0.1 N, 2.5 N \* 0.05 N.

**Technical Specification :**

Dimension (L\*W\*H CM): 16\*3.5\*2



**Laboratory Instrument India**