



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tiny Pressure Meter

**Product Code :**  
LBNY-0004-101000247



**Description :**

It is mainly used to demonstrate liquid and gas pressure in high school physics teaching. It is composed of a U-shaped tube, a scale plate, three through pipe, latex tubes, rod etc.

**Technical Specification :**

Dimension (L\*W\*H CM): 38\*4\*1.5



**Laboratory Instrument India**