



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Dual Cylinder

**Product Code :**  
LBNY-0004-101000252



**Description :**

The product for foot, hand press type dual-purpose, by suction nozzle, a cylinder cover, piston rod, a piston, a cylinder, a sealing ring, a bracket, an air inflating valve, outlet valve plate, pressure plate, a spring, a base and other components. The limit of the suction pressure is less than  $6.7 * 10^3 \text{Pa}$  (500mmHg), the minimum air pressure is over  $2.9 * 10^5 \text{Pa}$  (3kgf/cm<sup>2</sup>). Products by the outer cylinder and the suction nozzle, an air inflating valve, a base made of ABS engineering plastic pressing as a whole, the outer cylinder size diameter 36mm x 130mm, bottom diameter is 80mm.

**Technical Specification :**

Dual Cylinder



**Laboratory Instrument India**