



Email : sales@laboratoryinstrumentindia.com

Product Name :
Dual Cylinder

Product Code :
LBNY-0004-101000252



Description :

The product for foot, hand press type dual-purpose, by suction nozzle, a cylinder cover, piston rod, a piston, a cylinder, a sealing ring, a bracket, an air inflating valve, outlet valve plate, pressure plate, a spring, a base and other components. The limit of the suction pressure is less than $6.7 * 10^3 \text{Pa}$ (500mmHg), the minimum air pressure is over $2.9 * 10^5 \text{Pa}$ (3kgf/cm²). Products by the outer cylinder and the suction nozzle, an air inflating valve, a base made of ABS engineering plastic pressing as a whole, the outer cylinder size diameter 36mm x 130mm, bottom diameter is 80mm.

Technical Specification :

Dual Cylinder



Laboratory Instrument India