



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
High Frequency Signal Generator

**Product Code :**  
LBNY-0004-101000253



**Description :**

1. High-frequency signal frequency range: 400KHz -130MHz; (have 5 range)
2. Low-frequency signal: frequency 1KHz  $\pm$  5%, the rate of not less than 20mV continuously adjustable; power supply: AC 220V  $\pm$  10% 50Hz  $\pm$  5%.

**Technical Specification :**

High Frequency Signal Generator



**Laboratory Instrument India**