



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
The Experimental Apparatus Can Transform

**Product Code :**  
LBNY-0004-101000260



**Description :**

It is modular structure, and each part is placed in a special wooden box. This instrument is mainly used to do the experiments of the conversion between the mechanical energy, wind energy, electricity, light and heat.

**Technical Specification :**

Dimension (L\*W\*H CM): 13.5\*10.3\*6.8



**Laboratory Instrument India**