



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bedside Locker Psychiatric

Product Code :
LBNY-0001-820004



Description :

Bedside Locker Psychiatric

Technical Specification :

The sides and back of the cabinet of these bedside locker are made in one piece.
Top and bottom are welded to it. Door of the cabinet made of 1.25 mm thick mild steel sheet.
Cabinet door provided with a set of three ventilating louvers.
Legs fitted with rubber shoes. Pre treated and epoxy powder coated.
Unit a one piece seamless rotationally moulded design (for easy cleaning and disinfection)
Constructed from polyurethane or material of similar characteristics.
Features to allow base to be anchored to the floor (using bolts/screws/fasteners)
Base have 0mm clearance from ground all-round/no legs/no castors
Domed/rounded edges to ensure safety of patients and caregivers
Unit have fire retardant properties.
Dimension: D 400 x W 500 x H 800mm +/- 100mm.
Weight: ?15 kg.



Laboratory Instrument India

