



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mirror Psychiatric

Product Code :
LBNY-0001-820005



Description :

Mirror Psychiatric

Technical Specification :

The tamper proof PVC frame design is ligature free, resistant to breakage and easy to install. The mirror face is placed into the rebate from the back, then the ACM backing is glued into its own rebate. This method allows the mirror face to float slightly, however once fixed to a flat wall, the mirror and backing are totally encapsulated and cannot be removed without destruction of the frame (ie: the mirror cannot be pulled out from the front).

Wall Mounted Mirror With:

Constructed from highly polished stainless steel with Corian/or material of equivalent properties and safety features for psychiatry setup

Tamper proof design

Ligature resistant

Domed/rounded edges to ensure safety of patients and caregivers

Dimensions: Height 400mm x 600mm +/- 50mm.

Laboratory Instrument India

