



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Dining Room Table 6 Seater Capacity

**Product Code :**  
LBNY-0001-820006



**Description :**

Psychiatric Dining Room Table 6 Seater Capacity

**Technical Specification :**

This style of wooden dining table set creates a beautiful impression with the perception of modernization. The product requires basic self-assembly at the customer's end and comes with self-assembly instructions along with necessary accessories.

This product is a great addition to your living space which adds a touch of elegance to your modern decor theme.

Psychiatric Dining Room Table With:

One piece 6 seater dining table, seats and table top attached to the same frame.

Table top filled with high density dampening foam or equivalent material of similar properties.

Table top constructed from bi-layer and temper proof stainless steel surface. fastened with tamper resistant stainless steel screws OR equivalent design and materials.

The main frame constructed from heavy duty mild steel.

The frame capacity to be boiled to the floor.

ABS plastic seats or equivalent material (scratch resistant)

Dimensions:

Table top: (800mm +/- 50mm) x (2000mm +/- 200) X (35mm +/- 5mm)

Ground to table top height: 700mm +/- 15mm

Ground to seat height: 500mm +/- 50mm



**Laboratory Instrument India**