



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital Logic Trainer

**Product Code :**  
LBNY-0002-620011



**Description :**

Digital Logic Trainer

**Technical Specification :**

The Digital Trainer a modular course covering concepts, theory and applications for digital electronics training. It designed start the student with the basics of digital technology and goes onto advanced digital topics. All necessary equipment for digital trainer experiments such as multiple output power supply, clock generator, logic indicator, data switches, toggle switches, push buttons, pulsars and 7 segment displays installed on the base station.

The base station able to work with modules required for experiments.

A comprehensive operation manual and experiment manual are provided to be used together in doing the experiments.

Basic and Advance Modules

Fixed regulated DC power supplies : +5V, +15V, -15V, variable 1.5 ~ 15VDC

Logic indicators : 16 logic indicator

Data switches : 16 (for TTL & CMOS)

Toggle switches : 4 (SPDT)

Push buttons : 3 (2NO, 1NC)

Pulsers : 2 pulsers (A,A and B,B)

Seven segments : 4 seven segment digital displays

Clock generator : 1Hz ~ 100kHz

Test leads : 1 set of stackable 2mm test leads

---

Digital Logic Trainer, Digital Logic Trainer Equipments, Digital Logic Trainer Tools, Digital Logic Trainer Tool Kits, Digital Logic Trainer Manufacturers, Digital Logic Trainer Suppliers from India, China, Kenya



## Laboratory Instrument India