



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Resting and Seating Chair

**Product Code :**  
LBNY-0001-820008



**Description :**

Psychiatric Resting and Seating Chair

**Technical Specification :**

Recliner Chair has 6 Adjustable positions, Multi-Positions, All-weather to sit or sleep.

Cushion Cover zipper makes it easily removable for cleaning, Washable Cushion Cover.

Heavy-Duty Powder-coated steel frame and joints.

Load capacity up to 150kgs, No assembly required, Easy to Fold and Carry.

Unit a one piece seamless rotationally moulded design (for easy cleaning and disinfection) and avoid fluid ingress.

Constructed from Polyurethane or material of similar characteristics.

Ergonomically comfortable design

Ability feature to allow increase of unit weight (filled with sand/equivalent design). This feature prevents patients to use chair as a weapon.

Domed/ rounded edges to ensure safety of patients and caregivers

Four legged chair

unit have fire retardant properties.

Dimensions:

Height: 850mm +/- 25mm.

Width: 450mm +/- 25mm.

Weight: +/- 10kg.

Depth: 500+/- 50mm.



**Laboratory Instrument India**