




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Modular Ottomans

**Product Code :**  
LBNY-0001-820009



**Description :**

Psychiatric Modular Ottomans

**Technical Specification :**

Best Foot Rest – One-of-a-kind multi-purpose ottoman that allows you to keep your legs and feet up and relaxed for work from home long hours. Maintain a comfortable position while sitting at your desk or sofa, improves blood circulation, reduces swollen feet or ankles and helps in fast recovery from a leg injury.

Multipurpose – Can be used at home as a cube footrest stool, also can be used as a sitting stool for kids. The comfortable height perfect for resting your tired legs or feet as a foot rest stool, and as spare seat at the foot of bed, in living room or office.

Lightweight –Ottoman is made of High-quality Firm HR Foam which is light weight but Stable and easy to move around the house.

Ottomans with:

Rotationally moulded polyurethane base or material of similar characteristics.

Base mounted to the PVC top with tamper proof fasters or equivalent design.

Ergonomically comfortable design

Ability feature to allow increase of unit weight (filled with sand/equivalent design). This feature prevents patients to use chair as a weapon.

Interior/core made of flexible durable foam for- comfort and support/material of similar properties

Skin of the Ottoman made of flexible easy to clean and fire-retardant material (eg PVC).

Dimensions:

---

Height: 500mm+/- 50mm  
Width: 700mm +/- 50mm.  
Weight: +/- 10kg.



**Laboratory Instrument India**