



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sound Velocity Measurer

**Product Code :**  
LBNY-0004-101000321



**Description :**

It is a necessary instrument for measuring sound velocity in physical acoustics; it adopts advanced electronic circuit and high precision timer and can measure the speed of sound in the air at close range.

**Technical Specification :**

Dimension (L\*W\*H CM): 13\*8.5\*6.5



**Laboratory Instrument India**