



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sound Propagation Demonstrator

**Product Code :**  
LBNY-0004-101000324



**Description :**

It is mainly used for physics classroom demonstration experiment and student experiment in colleges and universities. It is used to demonstrate sound propagation in different mediums, such as air, liquid and solid, but not in vacuum.

**Technical Specification :**

Dimension (L\*W\*H CM): 24\*9\*34



**Laboratory Instrument India**