



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Wall Mounted Desk

**Product Code :**  
LBNY-0001-820010



**Description :**

Psychiatric Wall Mounted Desk

**Technical Specification :**

The strong sturdy table in Glossy White Colour can also be used as a marker/scribble board. Use it as a whiteboard.

Easy to clean due to the glossy surface. This desk can also be used as a standing desk.

Unit a one piece seamless rotationally moulded design (for easy cleaning and disinfection) and avoid fluid ingress.

Constructed from Polyurethane or material of similar characteristics.

Unit shelf space for extra storage.

feature to allow unit to be anchored to the wall using (Using bolts screws fasteners).

This feature prevents patients to use chair as a weapon.

Domed/ rounded edges to ensure safety of patients and caregivers

Unit have fire retardant properties.

**Dimensions:**

Height: 400mm+/- 50mm.

Width: 1000mm +/-100mm.

Depth: 400mm+/- 50mm.

Weight: ?10 kg.



**Laboratory Instrument India**