



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Psychiatric Modular Lounge Set

**Product Code :**  
LBNY-0001-820011



**Description :**

Psychiatric Modular Lounge Set

**Technical Specification :**

**Durable:** The Sofa Upholstery comes with non-woven fabric backing which becomes fuller and increases the tensile strength and dimensional stability. The sofa also has solid wood legs for greater strength and longer durability.

**Easy to move:** This lightweight sofa is easy to move around the house, so you can shift it to any room you wish without hassle. It also has good ground clearance, enabling you to easily clean underneath it.

**Stylish:** This sofa is lighter with a smooth finish, ensuring consistency and aesthetic good look that makes the furniture lightweight for easy portability. With high strength and stiffness to weight ratios, plywood is very cost-effective to use in structural applications.

Orthopaedic designed Lounge with multiple pieces collection with ability to reconfigure arrangement (L-shape, U-shape, Straight etc) based on therapy intentions to meet clients needs.

The configuration have corner pieces, seats with back rests and coffee table. Individual units have a weight bearing capacity of 250kg. Individual Pieces weight between 40kg and 50kg.

Minimum seating capacity of 8 clients. however, the lounge have ability to be extended to more capacity (pieces).

The core made of high density resilient, flexible, moulded polyurethane foam OR equivalent material properties. One piece moulded polyvinyl chloride cover (200-250mil thick or equivalent material).

Rounded edges and corners.

hall have high impact and resilient base. Units have floor anchoring/ganging straps accessories.

Piece have No springs, hard edges, hard corners, seams and padded fabric.

---

Dimensions:

Seat: height from ground: 450mm +/- 50mm, Depth: 600mm +/- 200mm.

Ottoman as Coffee Table: Height: 450mm +/- 50mm.

Depth: 250mm +/- 50mm.



**Laboratory Instrument India**