



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electricity Testing Set

**Product Code :**  
LBNY-0004-101000353



**Description :**

The pointer electroscope can test that glass rod is positively charged and the glue rod is negatively charged after the friction. The glass rod material is organic glass with silk, the glue rod material is polycarbonate with fur and the length is 300mm.

**Technical Specification :**

Dimension (L\*W\*H CM): 11.8\*7.4\*14.2



**Laboratory Instrument India**