



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ice Maker Cubes

Product Code :
LBNY-0001-830002



Description :

Ice Maker Cubes

Technical Specification :

Easy to clean construction with automatic cleaning function
User-friendly push button operation
Sturdy door improves efficiency and hygiene
Sturdy adjustable legs allow easy positioning
Manual defrost
Automatic filling - Professional water supply connection required
Soft bullet-shaped ice
Bench top
Output: 45-50kg/24hrs
Inbuilt storage bin: with cover
Storage bin capacity: 20-30kg
Automatic cut off when ice is full
Stainless steel finish
On/off switch
Makes ice cubes
With all plumbing accessories that makes it functional



Laboratory Instrument India