



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Ice Maker Cubes

**Product Code :**  
LBNY-0001-830002



**Description :**

Ice Maker Cubes

**Technical Specification :**

Easy to clean construction with automatic cleaning function  
User-friendly push button operation  
Sturdy door improves efficiency and hygiene  
Sturdy adjustable legs allow easy positioning  
Manual defrost  
Automatic filling - Professional water supply connection required  
Soft bullet-shaped ice  
Bench top  
Output: 45-50kg/24hrs  
Inbuilt storage bin: with cover  
Storage bin capacity: 20-30kg  
Automatic cut off when ice is full  
Stainless steel finish  
On/off switch  
Makes ice cubes  
With all plumbing accessories that makes it functional



**Laboratory Instrument India**