



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Demonstrating Adjustable Internal Resistance

**Product Code :**  
LBNY-0004-101000371



**Description :**

It is used for teaching demonstration in middle school physics; it can accurately do experiments, such as: closed circuit Ohm's law, using ammeter and voltmeter to measure electromotive force and internal resistance of battery, electromotive force and internal resistance of battery pack and output power of supply.

**Technical Specification :**

Dimension (L\*W\*H CM): 21.8\*8.4\*8.7



**Laboratory Instrument India**