



Email : sales@laboratoryinstrumentindia.com

Product Name :
Demonstrating Adjustable Internal Resistance

Product Code :
LBNY-0004-101000371



Description :

It is used for teaching demonstration in middle school physics; it can accurately do experiments, such as: closed circuit Ohm's law, using ammeter and voltmeter to measure electromotive force and internal resistance of battery, electromotive force and internal resistance of battery pack and output power of supply.

Technical Specification :

Dimension (L*W*H CM): 21.8*8.4*8.7



Laboratory Instrument India