



Email : sales@laboratoryinstrumentindia.com

Product Name :
Resistance Box

Product Code :
LBNY-0004-101000374



Description :

It is used to verify Ohm's law of part circuit; used as known resistance inserted in the circuit to regulate current; used as additional resistance to enlarge measurement range of ammeter; used to measure unknown resistance; used to measure electromotive force and internal resistance of power supply and is the essential instrument in electrical experiment.

Technical Specification :

Dimension (L*W*H CM): 16.6*13.2*7.8



Laboratory Instrument India