



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Upright Ultra-Low Freezer

Product Code :
LBNY-0001-830004



Description :

Laboratory Upright Ultra-Low Freezer

Technical Specification :

Features with microprocessor controlled system with interior digital temperature display to show various parameters simultaneously.

Designed with direct cooling system and manual defrost mode, equipped with eco- friendly refrigerant Mixture gas, Built-in high-precision microcomputer temperature control system and platinum resistor temperature sensors.

It offers high or low temperature alarm, Door open alarm and Keyboard lock, password protection function to prevent arbitrary adjustment of operating parameters.

Settable temperature in the range: -15 to -25°C

Temperature recorder

Forced air circulation

Audio and visual alarms

Double door, non-glass.

Capacity; 700 - 1000L

Internal material of construction – stainless steel

Shelves –6

Laboratory Instrument India

