



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laboratory Upright Refrigerator, Low Temperature Range

**Product Code :**  
LBNY-0001-830005



**Description :**

Laboratory Upright Refrigerator, Low Temperature Range

**Technical Specification :**

Seven temperature probes can ensure high precision of temperature control with nearly no fluctuation and thus can improve safety.

Equipped with a USB export interface, which can be used for storing data from the last month to the current month automatically in PDF format.

With a U-disk connected, the temperature data can be stored continuously and automatically.

The interior lighting system with double LED lights ensures high visibility inside the cabinet.

A test port is available for bringing convenience to users in testing temperature inside the cabinet.

Settable temperature in the range: 2 to 8°C

Temperature recorder

Audio and visual alarm for:

High/low temperature

Sensor and system failure

Power failure

Door open

Double door, glass.

Capacity; 700 - 1000L

Internal material of construction – stainless steel

Shelves – minimum. 12

Frostless design



**Laboratory Instrument India**