



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Demonstration Ohmmeter

**Product Code :**  
LBNY-0004-101000391



**Description :**

It consists of pointer-type internal magnetic structure and measuring circuit, with fourteen grades measurement range. It is used as galvanometer in student teaching demonstration experiment and is used to measure micro direct current, direct voltage and direct resistance.

**Technical Specification :**

Dimension (L\*W\*H CM): 27.1\*29.7\*11.8



**Laboratory Instrument India**