



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Demonstration Ammeter Voltmeter

**Product Code :**  
LBNY-0004-101000392



**Description :**

It consists of pointer-type internal magnetic structure and measuring circuit, with fourteen grades measurement range. It is used as galvanometer in teaching demonstration experiment and is used to measure direct current, direct voltage, AC current and AC voltage.

**Technical Specification :**

Dimension (L\*W\*H CM): 27.1\*29.7\*11.8



**Laboratory Instrument India**