



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Projection Current Detecting Meter

**Product Code :**  
LBNY-0004-101000394



**Description :**

It is used to measure voltage in a DC circuit in common middle school and technical secondary school teaching demonstration and experiment with projector. It can also be used as general voltmeter.

**Technical Specification :**

Dimension (L\*W\*H CM): 12\*10\*3.8



**Laboratory Instrument India**