



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mini Spin

Product Code :
LBNY-0001-830007



Description :

Mini Spin

Technical Specification :

Ideal for quick spin-downs and microfiltration
Swap rotor attachments can be interchanged without tools
Fast acceleration and braking
Exchangeable rotors
Tube Capacity: 6 x 1.5/2ml and 16 x 0.2ml PCR Strip
Maximum speed 7,000rpm, RCF: 2,680 xg
Run Time: Continuous
Driving Motor: DC Motor



Laboratory Instrument India