



Email : sales@laboratoryinstrumentindia.com

Product Name :
Projection Pneumatic Table

Product Code :
LBNY-0004-101000415



Description :

It is used to simulate Brownian motion, gas diffusion, particle scattering experiment, liquid evaporation, gas isothermal change, the cause of force changing object motion state, gas molecule in heat movement, motion of liquid molecular, uniform circular motion and gas compression-liquidation experiment.

Technical Specification :

Dimension (L*W*H CM): 29*29*6



Laboratory Instrument India